

## **MEDIUM FORMULATIONS**

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

### **53. COOKED MEAT MEDIUM**

Use commercially available Cooked meat medium (Oxoid CM 0081; BD226730) and supplement medium with

5 g/l yeast extract

5 g/l  $K_2HPO_4$

1 mg resazurin

Boil, cool, add 0,5 g/l cysteine and adjust pH to 7,0. Autoclave at 121°C for 15 minutes.