

MEDIUM FORMULATIONS

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

68. Photobacterium MPY MEDIUM

Peptone, Bacto (BD 211677)	5,0g
Yeast extract (BD 212750)	3,0g
NaCl	28,2g
MgSO ₄ . 7H ₂ O	6,9g
MgCl ₂ . 6H ₂ O	5,5g
CaCl ₂ . 2H ₂ O	1,5g
KCl	0,7g
Distilled water	1,0L

Adjust pH to 7,4. Autoclave at 121°C for 15 minutes.