

MEDIUM FORMULATIONS

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

36. Photobacterium BROTH

Tryptone	5,0 g
Yeast Extract	2,5 g
NaCl	30,0 g
NH ₄ Cl	0,3 g
MgSO ₄ . 7 H ₂ O	0,3 g
FeCl ₃	0,01 g
CaCO ₃	1,0 g
K ₂ HPO ₄	3,0 g
Na glycerophosphate	23,5 g
Agar	15,0 g
Distilled water	1000,0 ml

Adjust pH to 7,2. Autoclave 15 minutes at 121°C.