

## **MEDIUM FORMULATIONS**

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

### **55. PEPT-CARB-SOLUBLE STARCH**

Solution A:

Peptone	5,0g
Yeast extract	5,0g
K <sub>2</sub> HPO <sub>4</sub>	1,0g
MgSO <sub>4</sub> · 7H <sub>2</sub> O	0,2g
Soluble starch	20,0g
Agar	15,0g
Water	900,0ml

Solution B:

10% Na <sub>2</sub> CO <sub>3</sub> solution	100,0ml
--	---------

Autoclave the solutions separately at 115°C for 15 minutes, then mix Solutions A and B.