

MEDIUM FORMULATIONS

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

58. Escherichia MEDIUM

Trypticase peptone (BD 211921)	6,0g
K ₂ HPO ₄	0,2g
FeSO ₄ . 7H ₂ O	5,0mg
MgSO ₄ . 7H ₂ O	0,2g
L-Asparagine	0,15g
Glycerol	2,0ml
Vitamin B ₁₂	0,4mg
Agar (if needed)	15,0g
Distilled water	1,0L

Adjust medium for final pH 7,2. Autoclave at 121°C for 15 minutes.