

## **MEDIUM FORMULATIONS**

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

### **61. TYG MEDIUM**

Trypticase peptone (BD 211921)	10,0g
Yeast extract (Oxoid LP0021, BD 212750)	1,0g
Glucose	1,0g
NaCl	8,0g
CaCl <sub>2</sub> · 2H <sub>2</sub> O	0,3g
Distilled water	1,0L

Autoclave at 121°C for 15 minutes.