

## **MEDIUM FORMULATIONS**

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

### **63. LB MEDIUM WITH STREPTOMYCIN**

Tryptone (BD 211705)	10,0 g
Yeast extract (BD 212750)	5,0 g
NaCl	10,0 g
Streptomycin*	200,0mg
Distilled water	1,0L

Combine first three ingredients in 980 ml of water, adjust pH to 7,0 and autoclave at 121°C for 15 minutes.

Dissolve 200mg of streptomycin in 20 ml of water and to the cooled basal medium.