

## **MEDIUM FORMULATIONS**

The above media can be substituted by suitable analogues of variable origin. Instructions for the preparation of media are given as recommended by the Manufacturer.

### **64. FALCIVIBRIO MEDIUM**

Wilkins-Chalgren-Anaerobe-Broth (Oxoid CM0643) is supplemented with 0.3 g/l of L-cysteine and 1,0 mg/l of resazurin. Prepare the medium under 100% nitrogen gas. To the autoclaved medium add 50 ml/l of sterile serum (horse, bovine or lamb).